

## VM SERIES VIDEO MICROSCOPES "INFINITY CORRECTED OPTICAL SYSTEMS" F=200MM FOR INCIDENT (REFLECTED) LIGHT ILLUMINATION

| CAT NO.                                 | DESCRIPTION   | US\$ PRICE |
|---|---|------------|
| <b>VM SERIES VIDEO MICROSCOPES</b>      |   |            |
| VM-1V                                   | Video microscope, 1X magnification factor, with vertical "C" mount                | 865.00     |
| VM-1H                                   | Video microscope, 1X magnification factor, with horizontal "C" mount              | 1,035.00   |
| VM-1D                                   | Video microscope, 1X magnification factor, with vertical and horizontal "C" mount | 1,345.00   |
| VM-2V                                   | Video microscope, 0.5X magnification factor, with vertical "C" mount              | 895.00     |
| VM-035                                  | Video microscope, 0.35X magnification factor, with vertical "C" mount and iris    | 1,105.00   |
| <b>VM SERIES VIDEO MICROSCOPE STAND</b> |   |            |
| VM-PC                                   | Pole type stand w/ coarse and fine focus block for VM Series video microscopes    | 580.00     |

### **PLAN EPI METALLURGICAL OBJECTIVES FOR VM SERIES VIDEO MICROSCOPE**

| CAT NO.  | DESCRIPTION  | US\$ PRICE |
|--|--|------------|
| <b>BRIGHTFIELD PLAN EPI OBJECTIVES</b>                             |  |            |
| MA870  | Plan Epi Objective 5X, N.A. 0.10, W.D. 20.0 mm, F.N. 24      | 350.00     |
| MA871  | Plan Epi Objective 10X, N.A. 0.25, W.D. 7.48 mm, F.N. 24     | 450.00     |
| MA872  | Plan Epi Objective 20X, N.A. 0.40, W.D. 5.20 mm, F.N. 24     | 550.00     |
| MA873  | Plan Epi Objective 50X, N.A. 0.75, W.D. 0.38 mm, F.N. 24     | 650.00     |
| MA874  | Plan Epi Objective 100X, N.A. 0.90, W.D. 0.37mm, F.N. 24     | 1,150.00   |
| <b>BRIGHTFIELD LWD PLAN EPI OBJECTIVES (Long working distance)</b> |  |            |
| MA920  | LWD Plan Epi Objective 20X, N.A. 0.40, W.D. 11.0 mm, F.N. 24 | 850.00     |
| MA921  | LWD Plan Epi Objective 50X, N.A. 0.55, W.D. 7.5 mm, F.N. 24  | 1,050.00   |

### **ILLUMINATORS FOR USE WITH VM SERIES VIDEO MICROSCOPES**

| CAT NO.                                   | DESCRIPTION   | US\$ PRICE |
|---|---|------------|
| <b>150 WATT ILLUMINATOR FOR VM SERIES</b> |   |            |
| FT190/115                                 | Fiber optic light source, 115V AC, US cord, EKE lamp.   | 378.00     |
| FL150/115                                 | Fiber optic light source, 115V AC, US cord, EKE lamp.   | 378.00     |
| FL150/95                                  | Single arm light guide, flexible 36" for use with VM Series video microscope (7.9mm diameter at VM end, 15mm diameter at FL150/115 end) | 174.00     |
| FL150/70                                  | Halogen bulb, 21V 150W, for FL150/100 or FL150/200  | 38.00      |